Presentation:

Vial with 4 mL, 5 mL, 6 mL.

Indication:

It has the function of filling in the irregularities of the surfaces and allowing for the bond between the enamel and the restorative materials like, for example, compound resin and resinous cements.

Description of the Product:

The Power Bond 2.1 product incorporates two components, a primer and an adhesive agent, in a single vial, simplifying the clinical procedures, with excellent results in terms of adhesion forces. However, these easy-to-use and quick application systems require acid conditioning and multiple applications of adhesive to obtain a maximum adhesion to the enamel and the dentin.

Composition:

Bis-GMA, dimethacrylates, photoinitiator, ethanol, inert load, sodium fluoride.

Mode of use:

1. Perform the complete conditioning with phosphoric acid at 37% for 30 seconds on enamel and 15 seconds on dentin. Wash the area and then dry. Take care not to dehydrate the dentin when drying with air. The site should remain slightly moist.
2. Apply a fine layer of Power Bond 2.1 with the aid of a brush or burnisher. Apply a jet of air and polymerize for 20 seconds. Repeat the application if the surface is not smooth and shiny.
3. The site is prepared to receive the restorative material.
4. Perform the finishing of the restoration.

Precautions:

For exclusive use by dental surgeon.

Close the vial immediately after use.

Products such as eugenol or clove oil inhibit polymerization.

Do not use the product if the validity has expired.

Keep out of the reach of children.

Storage and Transportation Conditions:

Store the vial in a cool, dry location protected from dust, avoiding heat, sunlight, and other sources of direct light.

Do not store near materials that have eugenol or clove oil in their composition.

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