PRESENTATION: 01 syringe with 2g of NATURES VENEER cement; 05 application nozzles; Instruction manual for the professional.

COLORS: Yellow / Light / Clear / Yellow + / Light +

INDICATION: NATURES VENEER light-curing resin cement is indicated for the cementation of "dental contact lenses" and indirect facets that do not present opaque infrastructure. It provides permanent adhesion through the bonding and retention of the substrates through the photopolymerization system, being presented in a single component.

MODE OF USE:

1°To remove with brush/cup and prophylactic paste the temporary facet and the residues of temporary cement; 2°Select the color of the paste according to the color of the created facet and adjacent teeth; 3°Remove the part carefully and thoroughly wash the part and surface of the tooth with water.

Tooth preparation:

1°Check that the tooth surface and the prosthetic part do not contain any substance;
2°To make the relative or absolute isolation according to the professional preference;
3°Apply the acid conditioning of the enamel surface (30 s) and/or dentine (15 s) with CONDITIONING AGENT POWER ETCHING (BM4) or another of your choice. Then wash and dry without dehydrating;
4°If the cervical border is subgingival and if you prefer to use the relative insulation, do the procedure with the removal using retractor wires;
5°Apply two layers of dental adhesive rubbing with applicator for 10 seconds. Do not light-cure.

Cementation:

1°Do not handle the cement under direct light of the reflector, as this can accelerate the Natures Veneer polymerization; 2°Lay a thin layer of cement under the inner face of the facet; 3°Place the part in the preparation and, under slight pressure, remove the excess cement with exploratory probe or microapplicators; 4°Apply the oxygen inhibiting gel, such as POWER BLOCK (BM4) or another of your choice on the restoration margins after the removal of excess cement and before photoactivation, to avoid the layer inhibited by oxygen. Perform the photopolymerization for 40 seconds, on each of the margins (mesial, distal and cervical) and directly on the vestibular face; 5°Finish and polish with fine sandpaper discs, sandpaper strips, polishing pastes and felt strips.

* For photopolymerization, use halogen or LED light equipment with a power of 450 mW/cm² or more and a blue light emission spectrum (400-500 nm).

PRECAUTIONS: Exclusive use of the professional dentist; Follow the instructions carefully; Avoid contact with eugenol-based products as this may delay the product hardening; Keep out of the reach of children and domestic pets; Do not use the product if it is after the expiration date; Do not reuse empty packaging; Irritation may occur with the dental pulp, in case of contact with it or in places near it; In case of allergic reactions, discontinue use; If there is a pulpal exposure, protect with a protective material; Avoid contact with skin. In case of contact, remove excess immediately with water-absorbing cotton wool and rinse thoroughly with plenty of water; Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice; It is recommended to adjust the light of the reflector so that it does not directly affect the restoration, promoting a premature photopolymerization of the cement; Do not reuse the applicator nozzles to avoid cross-infection. Check correct syringe closure after use; In cases of allergic reactions to the product, discontinue use immediately.

STORAGE: Keep the product away from direct light and heat. Temperature: 5 - 30 °C.

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